

MEDICAL IRRIGATION TOWER

INTRODUCTION

OKLand[®] Medical Irrigation Tower can be applied in many application areas, such as orthopedics, obstetrics and gynecology, urology, emergency department and operating room. The nurse-friendly design can speed up the operation process and reduce the nurse's work intensity.

FEATURES

High Load Capability

- ☆ Totally up to 40,000 cc (10,000 cc max capacity per hook).
- ☆ Efficient continuous perfusion to avoid air entry into the joint cavity.

User-friendly

- ☆ Lift-Assist Technology: Easy to adjust the vertical distance.
- ☆ Four adjustable hooks: Can be operated independently.
- ☆ Hangers range: 1690 mm-2590 mm.
- ☆ Accurate Scale: Hanger height is displayed clearly both in cm and inch.
- ☆ Stainless steel handle: Hanging multiple injection pumps or other small medical devices, saving space.

Simple: Spring Design


- ☆ No hydraulic oil.
- ☆ No battery.

Stable

- ☆ Basement with universal locking wheel: Safe, stable and noiseless.

Economy


- ☆ Complete recyclability.
- ☆ No cost of consumables.





SPECIFICATIONS

Material:	Stainless steel
Minimum length/height (mm):	1,690
Maximum length/height (mm):	2,590
Method of height adjustment:	Lift-Assist Technology
Number of hooks:	4
Total load capacity (cc):	40,000
Load capacity per hook (cc):	10,000



Australian Distributor:



T: 1300 106 003

Dir: 07 3117 9539

E: sales@hmsmedical.com.au

www.hmsmedical.com.au